

LOT: 3 Immune affinity column SPE cartridge

No	Type of item	Unit	Quantity
1	Immuno affinity column, B1,G1,B2,G2, <ul style="list-style-type: none"> • Powerful: designed exclusively for laboratory use • Enhance recovery: improve aflatoxin G2 recovery • Exclusive: specifically for HPLC use • Comprehensive: determines individuals aflatoxin levels (B1, B2, G1 and G2) and total reading for all aflatoxins • Durable: long shelf life; requires no refrigeration • Versatile: for use with a variety of samples • Quick: 15 minutes to isolate toxin • Wide range: detects levels as high as 500ng and as low as 0.005 ng for total aflatoxin • Fast flow: passes more volume over the column • Capacity: 6 ml Pack of 25	packs	20
2	Supelclean LC-18 SPE tube bed weight500mg, volume 6ml, pack of 30	Pack	10
3	Supelclean LC-8 SPE tube bed weight500mg, volume 6ml, pack of 30	Pack	10
4	Supelclean LC-Ph SPE tube bed weight500mg, volume 6ml, pack of 30	Pack	10
5	Supelclean LC-CN SPE tube bed weight500mg, volume 6ml, pack of 30	Pack	10
6	Supelclean LC-SAX SPE tube bed weight500mg, volume 6ml, pack of 30	Pack	10
7	Supelclean LC- SCX SPE tube bed weight500mg, volume 6ml, pack of 30	Pack	10
8	Supelclean LC-WCX SPE tube bed weight500mg, volume 6ml, pack	Pack	10

	of 30		
9	Supelclean LC-NH2 SPE tube bed weight500mg, volume 6ml, pack of 30	Pack	10
10	Supelclean LC-Alumina SPE tube bed weight500mg, volume 6ml, pack of 30	Pack	10
11	Supelclean LC-Florisil SPE tube bed weight500mg, volume 6ml, pack of 30	Pack	10
12	Supelclean LC- Si SPE tube bed weight500mg, volume 6ml, pack of 30	Pack	10
13	Supelclean LC-diol SPE tube bed weight500mg, volume 6ml, pack of 30	Pack	10